# Dreamzon

## LightOn

*User Manual*

*Dreamzon Ltd.*

*E-mail: [info@dreamzon.com](mailto:info@dreamzon.com) [www.dreamzon.com](http://www.dreamzon.com)*

Copyright© 2010 DreamZon Ltd. All rights reserved.

Other product and company names mentioned herein may be trademarks or trade names of their respective owners. To the maximum extent permitted by applicable law, under no circumstances shall DreamZon or any of its licensors be responsible for personal injury or any incidental, special, indirect or consequential damages, howsoever caused. The third-party applications provided with your device may have been created and may be owned by persons or entities not affiliated with or related to DreamZon. DreamZon does not own the copyrights or intellectual property rights to the third-party applications. As such, DreamZon does not take any responsibility for end-user support, functionality of the applications, or the information in the applications or these materials. DreamZon does not provide any warranty for the third-party applications.

By using the applications you acknowledge that the applications are provided as is without warranty of any kind, express or implied, to the maximum extent permitted by applicable law.

Contact DreamZon at [info@dreamzon.com](mailto:info@dreamzon.com), or visit us at [www.dreamzon.com](http://www.dreamzon.com)

## **Safety Precautions and Important Notices**

- Do not try to connect cradle to electrical source! The LightOn product is operated by batteries and should not be connected to external DC power in any case.
- Do not use cradle on surface that vibrates (dryer, washing machine, etc...). LightOn should only be placed on flat non-vibrating surfaces, for example: table, shelf, countertop etc.
- Do not use near water – cradle is not waterproof. Keep the LightOn unit away from any heat source as it may damage the unit.

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

## About LightOn

LightOn is a unique decorative solution providing light notification for incoming SMS, Voice or Video calls, thus alerting its users that someone is trying to reach them. When a mobile phone is placed in the LightOn cradle, all incoming SMS messages or phone calls trigger a brightly flashing blue LED light. This light continues to flash until the device is reset, notifying users on the missed call.

The LightOn cradle's back-rest is designed to facilitate the use of sign language via mobile video phones. The angle of the cradle enables you to place your phone in the cradle and then comfortably sign with both hands.

LightOn allows simple connectivity to popular external signaling systems.

**Important Note:** Please check with your local Reseller the exact LightOn model, which supports your existing signaling system.

## **Table of Contents**

Safety Precautions and Important Notices	3
About LightOn	4
Getting Started / Unpacking	6
Device Layout	7
Battery Installation	9
Battery Level, Power Indicator	9
Using the LightOn Cradle	10
Preparing a Mobile Phone for Use	10
Light Panels	10
Increasing the LightOn Sensors Sensitivity	11
Battery Level Indicator	12
Resetting Light Indicator	12
Operating Environment	13
External Alerting and Triggering Systems	14
Troubleshooting	15

## Getting Started

### Unpacking

Your package contains the following items:

LightOn Cradle

Batteries

Cable

User Manual

Element	Function
Battery compartment	Holds the three AA batteries.
Service light Power Indicator	Indicates that LightOn cradle is working.
Reset button	Resets the light panel after a call or SMS has come in.
Light panel	Flashes repeatedly when the phone receives a call or SMS.
Back light panel	Flashes repeatedly when the phone receives a call or SMS.
Cable and connection socket	Connects LightOn to external triggering systems.

Service light Power Indicator

Light panel

Reset button

Cable and connection socket

Back light panel

Battery compartment

## Getting Started

### Battery Installation

1. Press firmly downwards on the battery case cover, on the sign. 
2. Insert three AA batteries, according to the (+) sign.
3. Click the battery case cover closed.

### Battery Level, Power Indicator

A blue LED light blinking on the top left of your LightOn cradle indicates how much charge is left in your batteries:

- An intermittent flashing every 15 seconds indicates that your batteries are properly charged.
- An incessant flashing indicates that battery power is weak. When battery power becomes too weak for the device to operate, it switches off automatically. Please replace the batteries if the Battery Level indicator flashes fast repeatedly, otherwise the LightOn device may not function well.

## Using the LightOn Cradle

### Preparing a Mobile Phone for Use

To use the LightOn cradle:

1. Place the mobile phone on vibrate mode. If a phone has multiple vibration modes, select the mode that uses the **strongest and longest vibration time**.
2. **Gently** place phone in cradle.

**Note:** If the phone is not placed gently, the movement may register as a vibration, causing the phone to light up. If this happens, press the Reset button.

### Light Panels

When the mobile phone on the LightOn cradle receives an SMS or phone call, a blue LED light begins flashing brightly in the front and back panels.

For the first two minutes after a call comes in, the light flashes rapidly. After two minutes, the frequency of the flashing decreases, conserving battery power, until it totally disappears.

## Increasing the Sensitivity of LightOn Sensors

The LightOn device works by converting the vibration signals of your cell phone into a flashing LED light. If your LightOn device is not able to recognize your cell phone's vibration signals (because the vibration is weak or short), it is possible to increase the sensitivity of LightOn's sensors.

### **How to Modify Sensor Sensitivity?**

1. Insert three AA batteries into the battery compartment. All the LED lights flash brightly, followed by the intermittent flashing of the left and right lights.
2. While the left and right LED lights are flashing, press the Reset button.
3. The LED lights will stop flashing and a single LED light will be lit on the upper right side of the cradle for a few seconds. LightOn is now properly calibrated and ready to be used.

**Note:** If you want to restore the default sensitivity settings of LightOn's sensors, remove one of the batteries from the battery compartment, hold the Reset button down for a few seconds and then replace the battery in its compartment.

**Note:** After changing the batteries, the LightOn device will return to its default sensitivity settings. If you want to restore the hi-sensitivity setting, repeat steps 1-3 above.

## **Battery Level Indicator**

The LED service light, on the upper left corner of the cradle, flashes every 15 seconds to indicate that the device is working. When battery power is low, this light flashes constantly, alerting you of the need to replace your batteries.

## **Resetting Light Indicator**

After receiving a call or SMS, you should reset the light indicator so that you know if a new call or SMS comes in.

### **To reset the light indicator**

- Click the Reset button in the front of the device. The device is now reset.

## **Operating Environment**

The LightOn phone cradle should be placed on a level, non-vibrating surface.

Since LightOn technology is based on vibration sensors, the cradle's operation is designed to have minimum disturbances of radio frequency signals or other mobile phones in the vicinity.


### **External Alerting and Triggering Systems**

LightOn can connect to external triggering and wireless signaling systems, thus enabling better communications and improving day to day quality of life.

LightOn supports today connectivity to market leading wireless signaling systems.

**To connect to an external alerting system:**

Connect one end of the provided cable into the LightOn cradle's connection socket, and the other end into the external system device connection socket. For example, connect LightOn to a transmitter, or a vibrating pillow.

## Troubleshooting

Problem	Possible Causes/Solutions
Light is not activated by SMS or phone call.	<p>Check that the batteries are inserted properly and that they do not need to be changed.</p> <p>Check that the phone is on the longest &amp; strongest vibrating mode.</p>
Light flashes even though no incoming call.	<p>Cradle placed on uneven surface.</p> <p>Cradle is placed very close to strong RF system.</p> <p>Push the reset button &amp; place the Mobile phone again</p> <p>Cradle placed on vibrating surface (for example, countertop with kitchen tools).</p> <p>Mobile phone was not placed gently in cradle – simply press the Reset button.</p>

Problem	Possible Causes/Solutions
Top left LED is constantly flashing.	Replace batteries.
I've connected the LightOn device to a wireless signaling system, but it doesn't trigger it.	Make sure that you've purchased the correct LightOn model . Please contact your reseller.

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Copyright© 2010 DreamZon Ltd. Alle Rechte vorbehalten.
Copyright© 2010 DreamZon Ltd. Tous droits réservés.
Copyright© 2010 DreamZon Ltd. Todos los derechos reservados.