

Vibrating Alarm Clock Model VA3 Quick Operating Guide

It's as easy as 1,2,3

1. To set time

- 1. Press and hold SET for 2 seconds
- 2. Press **UP** to display 12-hour clock or 24-hour clock; press **SET** to select
- Press UP or DOWN to set the Hour, press SET to select. Repeat these steps to set Minute.

2. To set alarm time & snooze

- Press MODE to go to Alarm mode, press and hold SET for 2 seconds to select
- Press UP or DOWN to set the alarm hour, press SET to select. Repeat these steps to set Minute.
- 3. Press UP to turn on or off the SNOOZE function, press SET to select

3. To set Alarm on/off & alarm mode

- Slide the switch to OFF to disable the alarm or to stop the on-going alarm action
- Slide the switch to AUDIO, SHAKER or ALL (audio + shaker) to select your alarm mode and to set the Alarm On

4. To review and test your alarm setting

1. Press and hold **DOWN** for 3 seconds, your alarm time appears and the alarm is triggered for a period of 5 seconds and then returns to the clock mode

5. To use clock as a Timer

- Press MODE to go to the Timer mode, press and hold SET for 2 seconds to select
- Press UP or DOWN to set the desired hour setting, press SET to select. Repeat these steps to set the desired minute setting
- 3. Press **UP** to start/stop the timer count-down
- When the set time expires, the alarm will sound-off (same alarm mode as the alarm clock setting) for 30 seconds (unless stopped by pressing any of the top keys)

6. To show the room temperature

- 1. Press **MODE** to go to the Thermometer mode, the room temperature appears
- 2. If needed, press UP to set reading in C or F
- 3. Press MODE to return to the clock mode

Need Help? Call Toll-free Help Line: (866)-376-9271 (9AM-5PM PST, M-F)